



DETERMINATION OF PESTICIDE RESIDUES AND CANNABIS POTENCY IN MARIJUANA AND CANNABIS FOODS

XIAOYAN WANG AND BRIAN KINSELLA

UCT, INC. 2731 BARTRAM ROAD, BRISTOL PA 19007, USA

CONTACT: XWANG@UNITEDCHEM.COM



BACKGROUND:

As of July 2015, 23 states and Washington D.C. in USA have legalized the medical use of marijuana, while 4 states and Washington D.C. have legalized the recreational use of marijuana. As a result, many testing labs are looking for fast, reliable, and cost-effective methods to determine cannabis potency and chemical residues, (e.g. pesticides, mycotoxins, and heavy metals) in marijuana and cannabis foods (more informally known as edibles). This poster utilizes the advantages of the QuEChERS (acronym for Quick, Easy, Cheap, Effective, Rugged, and Safe) to extract 35 pesticides and 3 cannabinoids, including tetrahydrocannabinol (THC), cannabidiol (CBD), and cannabiviol (CBN), in edibles and seized marijuana, followed by either serial dilutions for cannabis potency analysis, or a dispersive solid phase extraction (dSPE) cleanup for pesticide residue analysis.

The proprietary extraction salts are designed to enhance phase separation and partitioning of pesticides and cannabinoids into the extraction solvent layer, while the proprietary dSPE sorbent mixture removes the majority of matrix co-extractives including lipids, organic acids, sugars, artificial colors, and natural pigments (e.g. chlorophyll and anthocyanins) in various sample matrices tested in this study, resulting in clean extracts for instrumental detection.

EXPERIMENTAL:

QuEChERS Products	
Description	UCT Part Number
Pouches containing a proprietary blend of QuEChERS salts for THC Potency and Pesticide Testing; 50-mL centrifuge tubes included	ECQUUS10-MP
2-mL centrifuge tubes containing proprietary blend of dSPE sorbents for Pesticide Testing in Edibles and Marijuana	ECQUUS10-2CT

(a) Sample pre-treatment

- For hard candies, cookies, chocolate bars, and marijuana samples, grind to fine powders using a SPEX 6770 freezer mill.



Figure 1: Hard candy before (left) and after (right) freezer mill grinding

- For gummy samples, cut into slim pieces. Although freezer mill can grind gummies to powder at low temperature with the use of liquid nitrogen, it returns to gel state when temperature goes up to room temperature, thus gummy samples were cut instead of grinded.



Figure 2: Degassing of Reef cola (left) and Orange kush (right)

- For sodas, degas for 30 min by sonication.

(b) QuEChERS extraction

- Weigh 1 g of the pre-treated samples (hard candies, gummies, brownies, cookies, chocolate bars, oil, and marijuana) into 50-mL centrifuge tubes, add internal standard (optional) and 10 mL of reagent water, and hydrate for 1 hr using a horizontal shaker. For sodas, add 10 mL of the degassed sample and internal standard (optional) to 50-mL centrifuges.
- Add 10 mL of acetonitrile (MeCN) with 1% acetic acid.
- Add QuEChERS extraction salts from pouches (ECQUUS10-MP), and vortex for 10 sec to break up salt agglomerates.
- Shake for 1 min at 1000 stroke/min using a SPEX Geno/Grinder. For gummy samples, add 2 metal balls and shake for 10 min at 1000 stroke/min.
- Centrifuge at 3000 rcf for 5 min.




Figure 3: Samples after QuEChERS extraction (from left to right: hard candies, gummies, soda, chocolate, and marijuana)

(c) dSPE cleanup for pesticide residue analysis

- Transfer 1 mL of the supernatants to 2-mL dSPE tubes (ECQUUS10-2CT).
- Shake for 1 min at 1000 stroke/min using the SPEX Geno/Grinder.
- Centrifuge at 3000 rcf for 5 min.
- Transfer 200 μ L extract to the 2-mL auto-sampler vials, add 200 μ L of DI water, and vortex for 30 sec.



Figure 4: Comparison of QuEChERS extracts before and after dSPE cleanup (from left to right: hard candies, gummies, and marijuana)

(d) Make serial dilutions for cannabinoid analysis

- Perform serial dilutions (200 to 20,000 times depending on the cannabinoid concentration in different samples) of the QuEChERS extracts to 100 to 200 ppb.
- Spike the diluted samples with 50 and 150% of the target cannabinoids, which are used to quantify the cannabinoid concentration according to the standard addition method.

(e) Analyze by LC/MS/MS

- Analyze samples by LC/MS/MS (Thermo Scientific UltiMate 3000 LC system coupled to TSQ Vantage tandem MS) equipped with an UCT Aqueous C18 HPLC column (SLAQ100ID21-3UM).



RESULTS:

A. Pesticide residue analysis

Matrix-matched calibration curves were generated using post spiked blank extracts of homogenized green tea leaves to semi-quantify the pesticide residues in edibles and seized marijuana. Appropriate amounts of the pesticide spiking solutions (0.1 and 1 ppm) were added into the tea extracts to generate 6-point matrix-matched calibration curves with concentrations at 5, 10, 25, 50, 100, and 250 ng/mL (TPP as IS at 200 ng/mL). The responses were found to be linear ($R^2 > 0.99$) over the concentration range. The limit of quantitation (LOQ) of this method was found to be 5 ng/mL in the extract or 50 ng/g in the marijuana and edibles, and 5 ng/mL in the soda samples.

Table 1: Accuracy and precision of pesticides in spiked tea samples

Compound	Spiked at 10 ng/mL		Spiked at 50 ng/mL	
	Recovery, %	RPD (n=6)	Recovery, %	RPD (n=6)
Methamidophos	80	11	81	12
Azinphos	81	14	83	12
Acifluorfen	83	13	75	23
Oxydemeton-methyl	74	25	80	23
Dinoseb	95	15	75	14
Permethrin	57	26	59	30
Azinthion	105	25	87	12
Triethylphosphorothioate	97	14	82	14
Carbendazim	98	15	74	12
Spinosin	97	12	97	11
Fenoxypipron sulfone	171	11	108	15
Fenoxypipron sulfonide	99	14	96	16
Selenite	173	14	102	14
Carbaryl	91	20	103	14
Tebuconazole	105	9	105	17
Thiophan-methyl	76	7	78	8
Famphur	101	13	101	13
Flutriafol	92	14	95	10
Trifloxystrobin	105	11	99	12
Abamectin	99	24	95	13
DEET	105	30	87	12
Methidathion	102	23	115	14
Trifluralin	67	21	102	16
Bifenthrin	134	21	98	21
Permethrin	83	14	84	16
Acetamiprid	96	16	105	12
Sulfoxide	100	15	99	13
Triphenylethylene	85	2	87	5
Zenfenoxin	86	3	81	5
Quinazolin	92	4	92	3
Carbendazim	77	5	77	3
Permethrin	94	4	97	3
Ethion	92	3	92	5
Permethrin	87	8	88	6
Chlorpyrifos	90	9	93	9




Figure 5: Matrix-matched calibration curve of Diazinon ($R^2 = 0.9995$)

Table 2: Pesticide residues detected in edibles and seized marijuana

Sample	Detected pesticides (semi-quantitatively)
Orange (brownie)	14 ng/g Chlorpyrifos
Orange kush (soda)	10 ng/g Bifenthrin
Leaf tea (soda)	Not detected
CBD oil	1000 ng/g Methidathion
CBD oil	Not detected
Cookie and cream (bar)	Not detected
Cookie brownie	87 ng/g Bifenthrin
Mint milk chocolate	Not detected
Monkay bar	Not detected
Mint (hard candy)	Not detected
Mint (hard candy)	Not detected
Sour gummies	Not detected
Sour fruit (gummy)	Not detected
Marijuana sample 1	414 ng/g Methidathion
Marijuana sample 2	1000 ng/g DEET
Marijuana sample 3	100 ng/g DEET
Marijuana sample 4	120 ng/g DEET + 1386 ng/g Chlorpyrifos
Marijuana sample 5	172 ng/g Acetamiprid, 488 ng/g DEET, 6822 ng/g Chlorpyrifos
Marijuana sample 6	178 ng/g Carbaryl, 691 ng/g DEET, 71 ng/g Methidathion
Marijuana sample 7	

B. Cannabis potency determination by standard addition method

Example: Cookie and cream bar, labeled with 30 mg of THC in 45 grams (equals to 667 μ g/g) After QuEChERS extraction of 1 g of the ground cookie and cream sample into 10 mL MeCN, the concentration of THC in the supernatant will be 66.7 μ g/mL. Serial dilutions ($\times 10 \times 500$) were made to dilute the extract to about 133 ng/mL, then the diluted samples were spiked with 70 (about 50%) and 210 ng/mL (about 150%) cannabinoids. The peak areas were plotted against the diluted sample (0), 50% spiked (70 ng/mL) and 150% spiked (210 ng/mL) samples, a 3-point linear curve (Figure 6) was generated. The concentration in the diluted sample was calculated by dividing the intercept by the slope. With the calculated concentration, the peak areas were re-plotted (Figure 7) and a linear curve with R^2 of 0.9999 was obtained, indicating that the standard addition method is effective for accurate analyte quantitation.




Figure 6: Plot of peak area against the unspiked sample (0), and samples spiked at 50% (70 ng/mL) and 150% (210 ng/mL) of cannabinoids




Figure 7: Re-plot of peak area against the actual concentrations: unspiked sample (127 ng/mL), and samples spiked at 50% (197 ng/mL) and 150% (337 ng/mL) of cannabinoids

Calculations:

THC conc. in the diluted sample = $5501801/43372 = 127$ ng/mL
THC in the cookie and cream bar = 127 ng/mL $\times 500 \times 10$ mL/g $\times 45$ g = 29 mg (very close to the labeled 30 mg THC)

Table 3: Comparison of labeled and detected cannabinoids in edibles (unit: mg)

Edibles	CBD		CBN		THC	
	Labeled	Detected	Labeled	Detected	Labeled	Detected
CBD oil	500	493	<5	ND	6	12
Reef cola					10	7
Orange kush (soda)					10	6
Edible (brownie)					90	60
Nectarine hard candy					10	6
Mint (hard candy)					100	49
Cookie brownie	NA	ND	NA	ND	10	14
Cookie and cream					30	29
Monkey bar					100	69
Mint milk chocolate					100	74
Sour gummies					100	96
Sour fruit (gummy)					10	8
Sweet'n sour (gummy)	NA	28	NA	ND	100	31

Table 4: Cannabinoid (%) in seized marijuana samples

Marijuana	THC%	CBN%
Marijuana sample 1	1.1	0.4
Marijuana sample 2	0.9	0.3
Marijuana sample 3	0.9	0.2
Marijuana sample 4	2.3	0.2
Marijuana sample 5	1.5	0.3
Marijuana sample 6	0.9	0.2
Marijuana sample 7	1.1	0.4

CONCLUSIONS:

A fast and effective method was developed for the determination of pesticide residues and cannabis potency in edibles and seized marijuana. Pesticide residues and cannabinoids were extracted using a proprietary blend of QuEChERS salts followed by either a dSPE cleanup for pesticide analysis, or serial dilutions for cannabinoid potency test. Pesticides were found to be present in both edibles (sodas, oil, and brownie) and seized marijuana (6 of the 7 tested samples) at varied concentrations. The detected amounts of cannabinoids were compared to those listed on the labels of the cannabis infused food products, where a small portion of the tested products were accurately labeled, while others were either higher or lower than labeled.

ACKNOWLEDGEMENT

Keith Tucker (Vice President of Marketing at SPEX SamplePrep, LLC) is acknowledged for kindly providing the 6770 Freezer mill and 2010 Geno/grinder. Erik Swiatkowski (UCT) is thanked for his help in grinding samples using the SPEX 6770 freezer mill.