



APPLICATION NOTE


Lysing of Bacterial Cells in the Geno/Grinder

With kind permission of Russell H. Vreeland, West Chester University, January 2003

The Geno/Grinder was tested to determine if this technology could be used to lyse bacterial cells. Standard 96-well titer plates were used with 400-600 μ silica grinding beads (Molecular Biology Grade, cat. no. 2166). The delivery of the beads into each cell of the titer plate can be accomplished in a number of ways. In this case micropipette tips were filled to the mark with grinding beads, and each tip emptied into a titer plate well. This technique will deliver approximately 0.4 grams of silica beads per well. Other delivery systems are commercially available.

Experiments were performed with 250 – 500 ml of bacteria: in one case the gram-negative, salt-tolerant *Halomonas elongata*, and in the other, gram-positive *Bacillus*. Both organisms were grown to late log state at 37°C. with shaking, then harvested via centrifugation at 7000 RPM for 10 minutes. 0.4 grams of silica beads were added to each well, and the plate was sealed. The filled titer plates were then shaken in the Geno/Grinder at a setting of 1450 strokes per minute, in one-minute intervals for periods ranging from one minute to twelve minutes. Culture filtrate from the shaken plates was then measured for optical density at 260 nm. The average data for three separate "well sets" is shown in Figure 1.

Figure 1



∴ APPLICATION NOTE SP020:
Lysing of Bacterial Cells in the
Geno/Grinder

∴ APPARATUS:
Geno/Grinder®

∴ APPLICATION:
DNA/RNA and Other
Extractions



SPEX SamplePrep
65 Liberty St
Metuchen, NJ 08840 USA
Tel: 732-623-0465
Fax: 732-906-2492
E-mail: Sampleprep@spexcorp.com
www.spexsampleprep.com

European Headquarters
SPEX Europe
2 Dalston Gardens
Stanmore, HA7 1BQ, UK
Tel: +44 (0) 208 204 6656
Fax: +44 (0) 208 204 6654
E-mail: spexeurope@spex.com
Web: www.spexeurope.com

Conclusion

The Geno/Grinder is capable of lysing bacterial cells in 96-well titer plates with 400-600 μ silica beads as grinding media. With increased grinding time culture filtrate shows a definite increase in optical density, indicating release of nucleic acids during grinding. Grinding times of 6 to 9 minutes appear to be suitable for producing measurable amounts of nucleic acid.



Your Science is Our Passion.™

SPEX SamplePrep
65 Liberty St
Metuchen, NJ 08840 USA

Tel: 732-623-0465

Fax: 732-906-2492

E-mail: Sampleprep@spexcsp.com

www.spexsampleprep.com

European Headquarters

SPEX Europe

2-4 Dalston Gardens

Stanmore, HA7 1BQ, UK

Tel: +44 (0) 208 204 6656

Fax: +44 (0) 208 204 6654

E-mail: spexeurope@spex.com

Web: www.spexeurope.com